

Receptor-type tyrosine-protein phosphatase C

Catalog No: #AP86819



Package Size: #AP86819-1 50ug #AP86819-2 100ug #AP86819-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Receptor-type tyrosine-protein phosphatase C
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Leukocyte common antigen,L-CA,Lymphocyte antigen 5,Ly-5,T200,CD45,Ly-5
Accession No.	Uniprot ID:P06800
Uniprot	P06800
GeneID	19264;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ptprc

Note: This product is for in vitro research use only