

Profilin-2

Catalog No: #AP86823



Package Size: #AP86823-1 50ug #AP86823-2 100ug #AP86823-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Profilin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Profilin II
Accession No.	Uniprot ID:P35080
Uniprot	P35080
GeneID	5217;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PFN2

Note: This product is for in vitro research use only