

Phospholipase D1

Catalog No: #AP86824



Package Size: #AP86824-1 50ug #AP86824-2 100ug #AP86824-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Phospholipase D1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PLD 1,Choline phosphatase 1,Phosphatidylcholine-hydrolyzing phospholipase D1
Accession No.	Uniprot ID:Q9Z280
Uniprot	Q9Z280
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pld1

Note: This product is for in vitro research use only