

Peroxiredoxin-4

Catalog No: #AP86826



Package Size: #AP86826-1 50ug #AP86826-2 100ug #AP86826-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Peroxiredoxin-4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Antioxidant enzyme AOE372,Peroxiredoxin IV,Prx-IV,Thioredoxin peroxidase A0372,Thioredoxin-dependent peroxide reductase A0372,Rno2CysPrx04
Accession No.	Uniprot ID:Q9Z0V5
Uniprot	Q9Z0V5
GeneID	85274;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Prdx4

Note: This product is for in vitro research use only