

## Prostaglandin E synthase 2

Catalog No: #AP86827



Package Size: #AP86827-1 50ug #AP86827-2 100ug #AP86827-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Prostaglandin E synthase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Membrane-associated prostaglandin E synthase-2,mPGE synthase-2,Microsomal prostaglandin E synthase 2,mPGES-2,PGES2
Accession No.	Uniprot ID:Q66LN0
Uniprot	Q66LN0
GeneID	493639;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PTGES2

Note: This product is for in vitro research use only