

## Phosphatidylinositol-glycan-specific phospholipase D

Catalog No: #AP86830



Package Size: #AP86830-1 50ug #AP86830-2 100ug #AP86830-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Phosphatidylinositol-glycan-specific phospholipase D
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	PI-G PLD, Glycoprotein phospholipase D, Glycosyl-phosphatidylinositol-specific phospholipase D, GPI-PLD
Accession No.	Uniprot ID:Q8R2H5
Uniprot	Q8R2H5
GeneID	291132;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Gpld1

Note: This product is for in vitro research use only