

# Major prion protein homolog

Catalog No: #AP86834



Package Size: #AP86834-1 50ug #AP86834-2 100ug #AP86834-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Major prion protein homolog
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	65-21 protein,Acetylcholine receptor-inducing activity,ARIA,PR-LP,PRN-P
Accession No.	Uniprot ID:P27177
Uniprot	P27177
GeneID	396452;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PRNP

Note: This product is for in vitro research use only