

## Prostaglandin G/H synthase 2

Catalog No: #AP86836



Package Size: #AP86836-1 50ug #AP86836-2 100ug #AP86836-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Prostaglandin G/H synthase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Cyclooxygenase-2, COX-2, Mitogen-inducible PGHS, PHS II, Prostaglandin H2 synthase 2, PGH synthase 2, Prostaglandin-endoperoxide synthase 2, CEF-147
Accession No.	Uniprot ID:P27607
Uniprot	P27607
GeneID	396451;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: PTGS2

Note: This product is for in vitro research use only