

Prostaglandin E synthase 2

Catalog No: #AP86841



Package Size: #AP86841-1 50ug #AP86841-2 100ug #AP86841-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Prostaglandin E synthase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	GATE-binding factor 1,GBF-1,Microsomal prostaglandin E synthase 2,mPGES-2,Gbf1,Pges2
Accession No.	Uniprot ID:Q8BWM0
Uniprot	Q8BWM0
GeneID	96979;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ptges2

Note: This product is for in vitro research use only