

## Prostaglandin G/H synthase 1

Catalog No: #AP86842



Package Size: #AP86842-1 50ug #AP86842-2 100ug #AP86842-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Prostaglandin G/H synthase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cyclooxygenase-1,COX-1,Prostaglandin H2 synthase 1,PGH synthase 1,Prostaglandin-endoperoxide synthase 1,Cox-1,Cox1
Accession No.	Uniprot ID:P22437
Uniprot	P22437
GeneID	19224;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Ptgs1

Note: This product is for in vitro research use only