

Phosphoglycerate kinase 1

Catalog No: #AP86843

Package Size: #AP86843-1 50ug #AP86843-2 100ug #AP86843-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phosphoglycerate kinase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Pgk-1
Accession No.	Uniprot ID:P09411
Uniprot	P09411
GeneID	18655;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Pgk1

Note: This product is for in vitro research use only