

## POU domain, class 3, transcription factor 4

Catalog No: #AP86844



Package Size: #AP86844-1 50ug #AP86844-2 100ug #AP86844-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	POU domain, class 3, transcription factor 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Brain-specific homeobox/POU domain protein 4,Brain-4,Octamer-binding protein 9,Oct-9,Octamer-binding transcription factor 9,OTF-9,Brn-4,Brn4,Otf9
Accession No.	Uniprot ID:P62515
Uniprot	P62515
GeneID	18994;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Pou3f4

Note: This product is for in vitro research use only