

Retinoic acid receptor beta

Catalog No: #AP86847



Package Size: #AP86847-1 50ug #AP86847-2 100ug #AP86847-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Retinoic acid receptor beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	RAR-beta,Nuclear receptor subfamily 1 group B member 2,Nr1b2
Accession No.	Uniprot ID:P22605
Uniprot	P22605
GeneID	218772;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Rarb

Note: This product is for in vitro research use only