

## Peptidoglycan recognition protein 1

Catalog No: #AP86850



Package Size: #AP86850-1 50ug #AP86850-2 100ug #AP86850-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Peptidoglycan recognition protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Peptidoglycan recognition protein short,PGRP-S,PGLYRP,PGRP,TNFSF3L,SBBI68,UNQ639/PRO1269
Accession No.	Uniprot ID:O75594
Uniprot	O75594
GeneID	8993;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PGLYRP1

Note: This product is for in vitro research use only