

Prohibitin-2

Catalog No: #AP86851



Package Size: #AP86851-1 50ug #AP86851-2 100ug #AP86851-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Prohibitin-2 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Chicken |
| Other Names | RCJMB04_1i23 |
| Accession No. | Uniprot ID:Q5ZMN3 |
| Uniprot | Q5ZMN3 |
| GeneID | 771124; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PHB2

Note: This product is for in vitro research use only