

Phosphatidylinositol-binding clathrin assembly protein

Catalog No: #AP86852



Package Size: #AP86852-1 50ug #AP86852-2 100ug #AP86852-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phosphatidylinositol-binding clathrin assembly protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Clathrin assembly lymphoid myeloid leukemia,CALM,Calm,Fit1
Accession No.	Uniprot ID:Q7M6Y3
Uniprot	Q7M6Y3
GeneID	233489;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Picalm

Note: This product is for in vitro research use only