

Pentraxin-related protein PTX3

Catalog No: #AP86853



Package Size: #AP86853-1 50ug #AP86853-2 100ug #AP86853-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pentraxin-related protein PTX3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Pentaxin-related protein PTX3
Accession No.	Uniprot ID:Q0VCG9
Uniprot	Q0VCG9
GeneID	541148;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PTX3

Note: This product is for in vitro research use only