

## Poliovirus receptor-related protein 1

Catalog No: #AP86854



Package Size: #AP86854-1 50ug #AP86854-2 100ug #AP86854-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Poliovirus receptor-related protein 1  |
| Brief Description  | Recombinant Protein  |
| Host Species       | E.coli   |
| Purification       | Greater than 90% by SDS-PAGE   |
| Species Reactivity | Hu   |
| Other Names        | Herpes virus entry mediator C, Herpesvirus entry mediator C, HveC, Herpesvirus Ig-like receptor, HlgR, Nectin-1, HVEC, PRR1, CD111 |
| Accession No.      | Uniprot ID:Q15223  |
| Uniprot            | Q15223   |
| GeneID             | 5818;  |
| Tag Info           | His  |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)   |

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: PVRL1

Note: This product is for in vitro research use only