

Receptor-type tyrosine-protein phosphatase alpha

Catalog No: #AP86856



Package Size: #AP86856-1 50ug #AP86856-2 100ug #AP86856-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Receptor-type tyrosine-protein phosphatase alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Protein-tyrosine phosphatase alpha,Lrp
Accession No.	Uniprot ID:Q03348
Uniprot	Q03348
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ptpa

Note: This product is for in vitro research use only