

Pyruvate carboxylase, mitochondrial

Catalog No: #AP86857



Package Size: #AP86857-1 50ug #AP86857-2 100ug #AP86857-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pyruvate carboxylase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Pyruvic carboxylase,PCB
Accession No.	Uniprot ID:Q29RK2
Uniprot	Q29RK2
GeneID	338471;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PC

Note: This product is for in vitro research use only