

Pyruvate carboxylase, mitochondrial

Catalog No: #AP86858



Package Size: #AP86858-1 50ug #AP86858-2 100ug #AP86858-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pyruvate carboxylase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Pyruvic carboxylase,PCB
Accession No.	Uniprot ID:P11498
Uniprot	P11498
GeneID	5091;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PC

Note: This product is for in vitro research use only