

Placenta growth factor

Catalog No: #AP86859



Package Size: #AP86859-1 50ug #AP86859-2 100ug #AP86859-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Placenta growth factor |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Bovine |
| Other Names | PIGF,PLGF |
| Accession No. | Uniprot ID:Q9XS47 |
| Uniprot | Q9XS47 |
| GeneID | 280894; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PGF

Note: This product is for in vitro research use only