

Melanocyte protein PMEL

Catalog No: #AP86862



Package Size: #AP86862-1 50ug #AP86862-2 100ug #AP86862-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Melanocyte protein PMEL
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Melanocyte protein Pmel 17, Premelanosome protein, Retinal pigment epithelial-specific protein, Silver locus protein homolog, PMEL17, RPE1, SILV
Accession No.	Uniprot ID:Q06154
Uniprot	Q06154
GeneID	281487;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PMEL

Note: This product is for in vitro research use only