

POU domain, class 3, transcription factor 1

Catalog No: #AP86865



Package Size: #AP86865-1 50ug #AP86865-2 100ug #AP86865-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	POU domain, class 3, transcription factor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Octamer-binding protein 6, Oct-6, Octamer-binding transcription factor 6, OTF-6, POU domain transcription factor SCIP, Tst-1, Oct-6, Otf-6, Otf6, Scip
Accession No.	Uniprot ID:P20267
Uniprot	P20267
GeneID	192110;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Pou3f1

Note: This product is for in vitro research use only