

## POU domain, class 3, transcription factor 3

Catalog No: #AP86866



Package Size: #AP86866-1 50ug #AP86866-2 100ug #AP86866-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	POU domain, class 3, transcription factor 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Brain-specific homeobox/POU domain protein 1,Brain-1,Octamer-binding protein 8,Oct-8,Octamer-binding transcription factor 8,OTF-8,Brn-1,Brn1,Otf8
Accession No.	Uniprot ID:P31361
Uniprot	P31361
GeneID	18993;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Pou3f3

Note: This product is for in vitro research use only