

Proline-rich acidic protein 1

Catalog No: #AP86870



Package Size: #AP86870-1 50ug #AP86870-2 100ug #AP86870-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Proline-rich acidic protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Epididymis tissue protein Li 178,Uterine-specific proline-rich acidic protein,UPA,UNQ608/PRO1195
Accession No.	Uniprot ID:Q96NZ9
Uniprot	Q96NZ9
GeneID	118471;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PRAP1

Note: This product is for in vitro research use only