

Perilipin-2

Catalog No: #AP86872

Package Size: #AP86872-1 50ug #AP86872-2 100ug #AP86872-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Perilipin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Adipophilin, Adipose differentiation-related protein, ADRP, ADFP, ADRP
Accession No.	Uniprot ID:Q9TUM6
Uniprot	Q9TUM6
GeneID	280981;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PLIN2

Note: This product is for in vitro research use only