

Cytosolic phospholipase A2

Catalog No: #AP86880

Package Size: #AP86880-1 50ug #AP86880-2 100ug #AP86880-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cytosolic phospholipase A2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	cPLA2,Phospholipase A2 group IVA,CPLA2,PLA2G4
Accession No.	Uniprot ID:P47712
Uniprot	P47712
GeneID	5321;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PLA2G4A

Note: This product is for in vitro research use only