

Pappalysin-1

Catalog No: #AP86882



Package Size: #AP86882-1 50ug #AP86882-2 100ug #AP86882-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pappalysin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Insulin-like growth factor-dependent IGF-binding protein 4 protease,IGF-dependent IGFBP-4 protease,Pregnancy-associated plasma protein A,PAPP-A
Accession No.	Uniprot ID:Q8R4K8
Uniprot	Q8R4K8
GeneID	18491;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pappa

Note: This product is for in vitro research use only