

Retinoic acid receptor responder protein 2

Catalog No: #AP86885



Package Size: #AP86885-1 50ug #AP86885-2 100ug #AP86885-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Retinoic acid receptor responder protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Chemerin
Accession No.	Uniprot ID:Q29RS5
Uniprot	Q29RS5
GeneID	508990;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: RARRES2

Note: This product is for in vitro research use only