

Patatin-like phospholipase domain-containing protein 2

Catalog No: #AP86886

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP86886-1 50ug #AP86886-2 100ug #AP86886-3 1mg

Description

Product Name	Patatin-like phospholipase domain-containing protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Adipose triglyceride lipase,Atgl
Accession No.	Uniprot ID:P0C548
Uniprot	P0C548
GeneID	100911615;361676;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pnpla2

Note: This product is for in vitro research use only