

85/88 kDa calcium-independent phospholipase A2

Catalog No: #AP86887



Package Size: #AP86887-1 50ug #AP86887-2 100ug #AP86887-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	85/88 kDa calcium-independent phospholipase A2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Cal-PLA2,iPLA2,Group VI phospholipase A2,GVI PLA2
Accession No.	Uniprot ID:P97570
Uniprot	P97570
GeneID	360426;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pla2g6

Note: This product is for in vitro research use only