

Pyrroline-5-carboxylate reductase 1, mitochondrial

Catalog No: #AP86890



Package Size: #AP86890-1 50ug #AP86890-2 100ug #AP86890-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pyrroline-5-carboxylate reductase 1, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	P5C reductase 1
Accession No.	Uniprot ID:Q58DT4
Uniprot	Q58DT4
GeneID	539606;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PYCR1

Note: This product is for in vitro research use only