

Pleckstrin homology-like domain family A member 2

Catalog No: #AP86895



Package Size: #AP86895-1 50ug #AP86895-2 100ug #AP86895-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pleckstrin homology-like domain family A member 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Imprinted in placenta and liver protein,Protein 50C15,lpl,Tssc3
Accession No.	Uniprot ID:O08969
Uniprot	O08969
GeneID	22113;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Phlda2

Note: This product is for in vitro research use only