

## GTP-binding nuclear protein Ran

Catalog No: #AP86899



Package Size: #AP86899-1 50ug #AP86899-2 100ug #AP86899-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	GTP-binding nuclear protein Ran
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	GTPase Ran,Ras-like protein TC4,Ras-related nuclear protein
Accession No.	Uniprot ID:P62825
Uniprot	P62825
GeneID	100499482;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:RAN

Note: This product is for in vitro research use only