

# V(D)J recombination-activating protein 1

Catalog No: #AP86902



Package Size: #AP86902-1 50ug #AP86902-2 100ug #AP86902-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	V(D)J recombination-activating protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	RAG-1
Accession No.	Uniprot ID:P24271
Uniprot	P24271
GeneID	423164;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:RAG1

Note: This product is for in vitro research use only