

Dolichyl-diphosphooligosaccharide--protein glycosyltransferase 48 kDa subunit

Catalog No: #AP86904

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP86904-1 50ug #AP86904-2 100ug #AP86904-3 1mg

Description

Product Name	Dolichyl-diphosphooligosaccharide--protein glycosyltransferase 48 kDa subunit
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	DDOST 48 kDa subunit
Accession No.	Uniprot ID:A6QPY0
Uniprot	A6QPY0
GeneID	510682;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DDOST

Note: This product is for in vitro research use only