

Platelet-activating factor acetylhydrolase IB subunit gamma

Catalog No: #AP86906

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP86906-1 50ug #AP86906-2 100ug #AP86906-3 1mg

Description

Product Name	Platelet-activating factor acetylhydrolase IB subunit gamma
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	PAF acetylhydrolase 29 kDa subunit,PAF-AH 29 kDa subunit,PAF-AH subunit gamma,PAFAH subunit gamma,PAFAHG
Accession No.	Uniprot ID:Q29460
Uniprot	Q29460
GeneID	282515;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PAFAH1B3

Note: This product is for in vitro research use only