

Proliferation-associated protein 2G4

Catalog No: #AP86907



Package Size: #AP86907-1 50ug #AP86907-2 100ug #AP86907-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Proliferation-associated protein 2G4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	IRES-specific cellular trans-acting factor 45 kDa,ITAF45,Mpp1,Proliferation-associated protein 1,Protein p38-2G4,Ebp1,Plfap
Accession No.	Uniprot ID:P50580
Uniprot	P50580
GeneID	18813;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pa2g4

Note: This product is for in vitro research use only