

Peroxisome proliferator-activated receptor alpha

Catalog No: #AP86908



Package Size: #AP86908-1 50ug #AP86908-2 100ug #AP86908-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Peroxisome proliferator-activated receptor alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PPAR-alpha,Nuclear receptor subfamily 1 group C member 1,NR1C1,PPAR
Accession No.	Uniprot ID:Q07869
Uniprot	Q07869
GeneID	5465;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PPARA

Note: This product is for in vitro research use only