

Pulmonary surfactant-associated protein C

Catalog No: #AP86912



Package Size: #AP86912-1 50ug #AP86912-2 100ug #AP86912-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pulmonary surfactant-associated protein C
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	SP-C,Pulmonary surfactant-associated proteolipid SPL(Val),SP5,SFTP2
Accession No.	Uniprot ID:P11686
Uniprot	P11686
GeneID	6440;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SFTPC

Note: This product is for in vitro research use only