

Phosphotriesterase-related protein

Catalog No: #AP86913



Package Size: #AP86913-1 50ug #AP86913-2 100ug #AP86913-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phosphotriesterase-related protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Parathion hydrolase-related protein
Accession No.	Uniprot ID:Q60866
Uniprot	Q60866
GeneID	19212;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Pter

Note: This product is for in vitro research use only