

Prostaglandin-H2 D-isomerase

Catalog No: #AP86921



Package Size: #AP86921-1 50ug #AP86921-2 100ug #AP86921-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Prostaglandin-H2 D-isomerase |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | Beta-trace protein,Cerebrin-28,Glutathione-independent PGD synthase,Lipocalin-type prostaglandin-D synthase,Prostaglandin-D2 synthase,PGD2 synthase,PDS |
| Accession No. | Uniprot ID:P41222 |
| Uniprot | P41222 |
| GeneID | 5730; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PTGDS

Note: This product is for in vitro research use only