

Serine/threonine-protein kinase N3

Catalog No: #AP86931



Package Size: #AP86931-1 50ug #AP86931-2 100ug #AP86931-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serine/threonine-protein kinase N3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Protein kinase PKN-beta, Protein-kinase C-related kinase 3, PKNBETA
Accession No.	Uniprot ID:Q6P5Z2
Uniprot	Q6P5Z2
GeneID	29941;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PKN3

Note: This product is for in vitro research use only