

RING-box protein 2

Catalog No: #AP86937



Package Size: #AP86937-1 50ug #AP86937-2 100ug #AP86937-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RING-box protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Rbx2,CKII beta-binding protein 1,CKBBP1,RING finger protein 7,Regulator of cullins 2,Sensitive to apoptosis gene protein,RBX2,ROC2,SAG
Accession No.	Uniprot ID:Q9UBF6
Uniprot	Q9UBF6
GeneID	9616;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:RNF7

Note: This product is for in vitro research use only