

Receptor-binding cancer antigen expressed on SiSo cells

Catalog No: #AP86938

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP86938-1 50ug #AP86938-2 100ug #AP86938-3 1mg

Description

Product Name	Receptor-binding cancer antigen expressed on SiSo cells
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cancer-associated surface antigen RCAS1,Estrogen receptor-binding fragment-associated gene 9 protein,RCAS1
Accession No.	Uniprot ID:O00559
Uniprot	O00559
GeneID	9166;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:EBAG9

Note: This product is for in vitro research use only