

1-acyl-sn-glycerol-3-phosphate acyltransferase alpha

Catalog No: #AP86939



Package Size: #AP86939-1 50ug #AP86939-2 100ug #AP86939-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	1-acyl-sn-glycerol-3-phosphate acyltransferase alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	1-acylglycerol-3-phosphate O-acyltransferase 1,1-AGP acyltransferase 1,Lysophosphatidic acid acyltransferase alpha,LPAAT-alpha,Protein G15,G15
Accession No.	Uniprot ID:Q99943
Uniprot	Q99943
GeneID	10554;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:AGPAT1

Note: This product is for in vitro research use only