

1-acyl-sn-glycerol-3-phosphate acyltransferase alpha

Catalog No: #AP86940



Package Size: #AP86940-1 50ug #AP86940-2 100ug #AP86940-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	1-acyl-sn-glycerol-3-phosphate acyltransferase alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	1-acylglycerol-3-phosphate O-acyltransferase 1,1-AGP acyltransferase 1,Lysophosphatidic acid acyltransferase alpha,LPAAT-alpha
Accession No.	Uniprot ID:O35083
Uniprot	O35083
GeneID	55979;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Agpat1

Note: This product is for in vitro research use only