

Parkin coregulated gene protein

Catalog No: #AP86945



Package Size: #AP86945-1 50ug #AP86945-2 100ug #AP86945-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Parkin coregulated gene protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Molecular chaperone/chaperonin-binding protein,PARK2 coregulated gene protein,GLUP
Accession No.	Uniprot ID:Q96M98
Uniprot	Q96M98
GeneID	135138;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PACRG

Note: This product is for in vitro research use only