

Progesterin and adipoQ receptor family member 3

Catalog No: #AP86947



Package Size: #AP86947-1 50ug #AP86947-2 100ug #AP86947-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Progesterin and adipoQ receptor family member 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Progesterin and adipoQ receptor family member III,Raf kinase trapping to Golgi,RKTG
Accession No.	Uniprot ID:Q6TCH7
Uniprot	Q6TCH7
GeneID	152559;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PAQR3

Note: This product is for in vitro research use only